Simply the best fixation.

Medical articulated arm with quick fixation + components
Medizin-Gelenkarm mit Schnellfixierung + Komponenten
Bras articulé médical avec fixation rapide + composants
Braccio articolato con fissaggio rapido per uso medicale + componenti

Instructions for use for non-sterilizable products
Gebrauchsanweisung für nicht sterilisierbare Produkte
Mode d'emploi pour les produits non stérilisables
Istruzioni per l'uso per prodotti non sterilizzabili

English

Introduction 4
 General Information..... 4
 Symbols used..... 4
Intended Use 4
Product description 5
Installation and use..... 5
 First check 5
 Warranty..... 5
 Handling of the base component..... 5
 Handling of the articulated arm..... 6
 Handling of the head components..... 7
Processing (cleaning and disinfection)..... 7
 Manual cleaning and disinfection by wiping – standard procedure 8
 Manual cleaning and disinfection - method for Art. 5.427 & 5.428 8
Control 9
Packaging and storage..... 9
Maintenance 9
Material resistance..... 9
Reusability..... 10

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Introduction

Thank you for buying our product. You have chosen a high quality Swiss product. We are confident that this product will meet your expectations.

General Information

Please read these instructions carefully. Improper handling can cause harm to the patient or damage the products. Use universal precautions when handling contaminated or biohazardous components / materials.

Manufacturer

Baitella AG
Thurgauerstr. 70
CH-8050 Zurich
Switzerland

Phone: +41 (0)44 305 80 00
Fax: +41 (0)44 305 80 05
E-Mail: info@fisso.com
Internet: www.fisso.com

Copyrights

All rights reserved. Duplication, adaptation, or translation, of any part of this document, without previous written authorization by Baitella AG is prohibited, except within the framework of the copyright regulations.

Reserved rights to technical changes

The drawings in this instruction may differ slightly from actual products.

Symbols used

Symbol	Definition	Danger	Consequence
	DANGER!	Immediate danger for people	Death or serious injuries
	WARNING!	Possible danger for people or objects	Health damages or severe material damages
	ATTENTION!	Possible danger to objects	Material damages

Symbol	Description
	Symbol for products that are put into use according to the directive 93/42 EEC
	NOTE Additional hints or important information
	See instructions for use
	Manufacturer

Intended Use

This product is used to hold and position breathing circuits, pressure transducers plates and other devices. It must only be used by trained medical personnel capable of judging and controlling any danger to patients. If not the case, the user assumes all responsibility.

DANGER!

The product will be delivered non sterile and is **non-sterilizable**. Prior to the initial and each following use, the product has to be cleaned and disinfected as well as checked for visual irregularities and malfunctioning according to the indications given in this manual.

Guideline for medical devices (conformity declaration)
This product corresponds to the guideline 93/42/EEC.

WARNING!

Do not attach height sensitive devices or transducers with drainage capacity. Injury to patient may occur.

For a „non-grounded“ patient positioning (complying to the safety norms EN IEC 60601-1 and 2-2, VDE 0750 Part 2 and VDE 0753) an insulated rail clamping base (Art. no. 5.240) must be used. Special care should be taken when these products are used together with high-frequency applications. Contacts between the product and the high-frequency device must be avoided.

Product description

The products are usually delivered as set.

A complete set includes:

- a clamping base (e.g. Rail clamping base or Universal clamp)
- a column (I, L, or U-shape)
- an articulated arm
- a head component (e.g. breathing circuit holder plate, clamp for pressure transducer plates)

The available standard sets and components can be viewed in the product brochures on the website (www.fisso.com) and can be downloaded. Individual FISSO sets can be assembled on customer demand.

Before delivery to the clients, the column and head components are normally rigidly threaded to the articulated arm and additionally fastened with glue or pins.

Installation and use

First check

Check the product immediately after receipt for eventual transport damages and for completeness. Complaints can only be considered if the manufacturer, distributor or carrier is immediately notified through damage protocol.

Warranty

The manufacturer assumes warranty for materials and manufacturing defects as well as proper functioning of FISSO holding systems for two years from date of purchase.

Not covered by the warranty are normal wear, the consequences of inappropriate handling, combinations with external products, or damage caused by the user. Other or further claims, in particular any claims for damage are excluded. The warranty applies only if the Baitella AG repair service is used. Please see also the General Delivery Conditions.

Handling of the base component

The set can be easily mounted to the surgical table or patient bed without additional tools. All handling should be done with as little force as possible and as much force as necessary!

	1	Attach the base to the rail of the surgical table or patient bed and secure it by turning the handle clockwise (Pos. 1) from below. Pay attention that the upper hook (Pos. 4) grabs behind the rail of the surgical table. Check if the base is fixed securely and firmly.
	Step 2 and 3 refer to bases with height adjustable columns	
	2	If necessary remove the stop lid (Pos. 3) from the column by turning counter clockwise. Place the column with one hand into the column opening of the rail clamping base and secure the column with the other hand by turning the handle (Pos. 2) clockwise. Then screw the stop lid into the column clockwise and tighten. Check if the column is fixed securely and firmly into the base.
3	The desired height and orientation of the column including the articulated arm can be adjusted at any time by loosening the handle (Pos. 2). Important: Hold the column with one hand and loosen the handle with the other hand. When the desired position is reached, tighten the handle.	

DANGER!


Risk of injury and damage to property
If the rail clamping base / column is not correctly tightened, it may become loose and may cause injuries to the patient.
The rail clamping base and if applicable column must be tightened properly.

ATTENTION!

The handle (Pos. 2) should not be tightened at any time without the column first being in place.

Handling of the articulated arm

	1	Hold the articulated arm with one hand in the anterior segment (A) and manipulate the central clamping handle (Z) with the other hand
	2	To release turn the central handle (Z) counter clockwise (B) as far as necessary
	3	Move the articulated arm into the desired position
	4	To secure , turn the central handle (Z) clockwise (C)
	5	Check that the articulated arm is fixed firmly and is functioning properly

DANGER!

Danger of serious injuries
The articulated arm can move unintentionally if the central clamping is released. Always hold the head component on the anterior segment (A) with one hand and manipulate the central clamping device (Z) with the other hand.

DANGER!

Danger of serious injuries
If the articulated arm is not tightened properly it can become loose and move, potentially causing injuries.
The articulated arm must be tightened properly.

DANGER!

Danger of serious injuries
A damaged product can cause injuries.
Use only products in proper condition and check functionality.

WARNING!

Risk of injury
The articulated arm can transmit electrical current and heat to the patient.
Avoid the contact between the articulated arm and any sources of electrical current or heat.

ATTENTION!

The fixation of the articulated arm is based on the principle of friction. Changing the position without loosening the clamping mechanism can cause damage and will shorten the lifespan of the articulated arm. The articulated arm can be adjusted with little force. If the central clamping handle (Z) has been fully loosened it has to be turned in the clockwise direction!

ATTENTION!

If the articulated arm is fixed to an apparatus, the apparatus must not be moved / shifted by means of the articulated arm!
Do not use the articulated arm to move / shift an apparatus.

Handling of the head components

The device (e.g. breathing circuit, pressure transducer) to be held by the head component must be secured safely.

DANGER!

Danger of serious injuries
If the device is not fixed properly it can come loose and move, which can cause injuries.
Check secure fixation of the device.

	<p>Handling of swinging head By opening and closing of the wing-screw (F) of the swinging head, the movability of the tube holder can be increased or decreased.</p>
---	---

Processing (cleaning and disinfection)

All products must be cleaned and disinfected prior to each application; this is also required for initial use after delivery (cleaning and disinfection after removal of the protective packaging). Effective cleaning and disinfection is an indispensable requirement for a safe application of the products.

You are responsible for the cleanliness of the products. Please pay attention to avoid severe contamination of the product during application (especially by handling with contaminated gloves) and follow the procedures below.

Additionally, please observe the legal provisions valid for your country as well as the hygienic instructions of the doctor's practice or of the hospital. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for the US).

Cleaning and disinfection

Depending on the product, select one of following cleaning and disinfection procedures:

- Manual cleaning and disinfection by wiping – standard procedure (for all unsterile applied products with the exception of single and twin tube holders)

- Manual cleaning and disinfection by dipping – method for plastic products with strongly profiled geometry (for single and twin tube holders only (Art.# 5.427 & 5.428), not permitted for all other products)

Pre-treatment

Please remove coarse impurities of the products directly after use (within a maximum of 2 hrs). To this end, wipe all visually contaminated surfaces of the products intensively with disinfection wipes¹.

¹ Please consider that the wipes should be alcohol- and aldehyde-free (otherwise there is a risk of fixation of blood impurities or damage to the product), possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), be suitable for disinfecting instruments made of metallic or plastic material and compatible with the products (see chapter "Material resistance"). Please consider that a disinfectant used in the pre-treatment step is solely for the safety of the staff and cannot replace the disinfection step later to be performed after cleaning.

Manual cleaning and disinfection by wiping – standard procedure

Observe the following points when selecting the wipes for cleaning and disinfection:

- Fundamental suitability for the cleaning and disinfection of instruments made of metallic or plastic material
- Approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking)
- Compatible with the products (see chapter Material resistance)

Pay attention to the instructions of the wipe manufacturer regarding application and contact time.

- Procedure:
1. Dismantle the system inasmuch possible.
 2. Wipe all surfaces of the product completely and carefully with fresh wipes. Visibly dirty or dry wipes must not be used further.
 3. Check all surfaces for remaining visible contamination and wipe again by applying a fresh wipe, if required.
 4. Check the products (see chapters "Control" and "Maintenance").
 5. Reassemble the product and wipe again using a fresh wipe.

Manual cleaning and disinfection - method for Art. 5.427 & 5.428

Observe the following points when selecting the cleaning and disinfection detergents:

- Fundamental suitability for the cleaning and disinfection of instruments made of plastic material
- Disinfectant with approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) compatible with the cleaning detergent used
- Compatibility of the used detergents with the products (see chapter "Material resistance")

Combined cleaning/disinfection detergents are not recommended.

Observe the instructions of the detergent manufacturers regarding concentration, temperature and soaking time as well as post-rinsing. Please use only freshly prepared solutions as well as only sterile or fully demineralised water and a soft, clean, lint-free cloth and/or filtered air for drying.

- Procedure:
- Cleaning
1. Remove the tube holder from the holding system.
 2. Soak the tube holder for the given soaking time in the cleaning solution so that the tube holder is sufficiently covered. Assist by careful brushing with a soft brush (never use metal brushes or steel wool), ultrasonic treatment is not recommended. Ensure that there is no contact with other products.
 3. Then remove the tube holder from the cleaning solution and post-rinse it at least three times with water.
 4. Check the tube holder (see chapters "Control" and "Maintenance").
- Disinfection
5. Soak the tube holder for the given soaking time in the disinfectant solution so that the tube holder is sufficiently covered. Ensure that there is no contact with other products.
 6. Then remove the tube holder from the disinfectant solution and post-rinse it at least five times with water.
 7. Dry and pack the tube holder immediately after removal (see chapter "Packaging")

DANGER!

The product will be delivered non sterile and is **non-sterilizable**. Prior to the initial and each following use, the product has to be cleaned and disinfected as well as checked for visual irregularities and malfunctioning according to the indications given in this manual.

ATTENTION!

Material damage due to improper cleaning and disinfection.
This product can be cleaned and disinfected **only manually**. The cleaning and disinfection agents must be compatible with the product.

WARNING!

The articulated arm **must not be immersed in liquids**.
The central handle of the articulated arm must be tightened during cleaning/disinfection.

Control

Check after cleaning/disinfection the action of moving parts to ensure smooth operation throughout the intended range of motion, the clamping force, the surface regarding corrosion, damages, chipping, and dirt and sort out impaired products. Products that are still dirty must be cleaned and disinfected again.

Packaging and storage

Please avoid re-contamination of the product by adequate protection / packaging and store the products in a dry and dust-free place. After assembly and prior to next use, further disinfection by wiping all surfaces (as in pre-treatment) may be required.

Maintenance

Special maintenance is not required. The product should be sent back to the manufacturer or distributor if the action of moving parts or the clamping force is impaired, or in case of damage.

Insulated rail clamping bases must not be used after the indicated expiry date but can be sent back to Baitella AG for service.

Instrument oil may never be used.

NOTE

Repair must be carried out by Baitella AG. If not respected, the warranty becomes void!

ATTENTION!

Please send back only reprocessed products to the distributor / manufacturer for repair (cleaned and disinfected).

Material resistance

Please note that the substances listed are **not** ingredients of the cleaning or disinfection agent:

- Organic, inorganic and oxidizing acids (permitted pH higher than 6.5)
- Bases (permitted pH lower than 8.5, neutral detergents on the basis of quaternary compounds or guanidine derivatives recommended)
- Organic compounds / solvents (for example: alcohol, aldehyde, acetone, ether, benzene)
- Oxidizing agents (for example: peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic, halogenated hydrocarbons

The product is non-sterilizable. Please do not expose any products to temperatures higher than 60 °C (140 °F).

Reusability

The products can be reused if processed properly and undamaged (see chapter "Control"). The user assumes full responsibility if damaged or dirty products are used (no liability will be accepted if this is not observed).

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