Simply the best fixation.

EN

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English

INTRODUCTION.....	16
GENERAL INFORMATION.....	16
SYMBOLS USED	16
INTENDED USE	18
PRODUCT DESCRIPTION	18
INSTALLATION AND USE.....	19
FIRST CHECK.....	19
WARRANTY.....	19
ASSEMBLY	19
DISASSEMBLY	20
HANDLING OF THE ARTICULATED ARM	20
HANDLING OF ARM SUPPORTS.....	21
PROCESSING (CLEANING AND DISINFECTION)	21
PRE-TREATMENT	21
MANUAL CLEANING AND DISINFECTION BY WIPING.....	22
CONTROL	23
PACKAGING AND STORAGE	23
MAINTENANCE.....	23
MATERIAL RESISTANCE	23
REUSABILITY.....	24
OBLIGATION TO NOTIFY	24
DISPOSAL	24
TECHNICAL DATA	24

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Introduction

Thank you for buying our product.
 You have chosen a high quality Swiss product.

General Information

Please read this instruction for use carefully and keep it close to the product. Improper handling can cause harm to the patient or damage the products. Use universal precautions when handling contaminated or biohazardous components / materials.

Manufacturer

Baitella AG
 Thurgauerstr. 70
 CH-8050 Zürich
 Switzerland

Phone: +41 44 305 80 00
 Fax: +41 44 305 80 05
 E-Mail: info@baitella.com
 Internet: www.fisso.com

Baitella GmbH
 Moltkestrasse 41
 DE-78532 Tuttlingen
 Germany

Copyrights

All rights reserved. Duplication, adaptation, or translation, of any part of this document, without previous written authorization by Baitella AG is prohibited, except within the framework of the copyright regulations.

Reserved rights to changes

Baitella AG reserves the right to make changes without notice in design, specifications and models and assumes no liability for the correctness and completeness of the content of the Instructions for Use.

Symbols used

Symbol	Definition	Danger	Consequence
	DANGER!	Immediate danger for people	Death or serious injuries
	WARNING!	Possible danger for people or objects	Health damages or severe material damages
	ATTENTION!	Possible danger to objects	Material damages

Symbol	Description
	Symbol for products that are put into use according to the regulation (EU) 2017/745 on medical devices
	NOTE Additional support or useful information
	See instructions for use
	Manufacturer
	Maximal carrying capacity of the articulated arm during intended use
	The articulated arm must not be sprayed directly with liquid or immersed in liquids. The central handle of the articulated arm must be tightened during cleaning / disinfection.
	European authorized representative
	Distributor
	Importer
	Medical device (in this instruction for use indicated as «product»)
	Unique Device Identification
	Article number
	Date of manufacture
	Lot number
	Serial number
	Quantity
	Gentle washing
	Do not bleach
	Iron with low heat
	Dry clean, any solvent except trichloroethylene

Intended Use

The Arm Support is used to position the patient arm during anaesthesia as well as during the initiation and reversal phase of anaesthesia. It must only be used by trained medical personnel capable of judging and controlling any danger to patients. If not the case, the user assumes all responsibility.

DANGER!

Unauthorized changes or modifications to the FISSO holding system are prohibited for safety reasons.

DANGER!

All FISSO components are designed for optimal performance when used together with the FISSO holding system. If a products / component from another manufacturer is used, the user takes full responsibility.

DANGER!

The product will be delivered non sterile and is **non-sterilisable**. Prior to the initial and each following use, the product has to be cleaned and disinfected as well as checked for visual irregularities and malfunctioning according to the indications given in this manual.

Guideline for medical devices (conformity declaration)

This product corresponds to the regulation (EU) 2017/745 on medical devices.

Special care should be taken when these products are used together with high-frequency applications. Contacts between the product and the high-frequency device must be avoided.

Product description

The product range along with all related information can be reviewed on www.fisso.com.

	Arm support lateral Art. # 6646.XX	Arm support lateral Art. # 4646.XX	Arm support vertical Art. # 6656.XX	Arm support vertical Art. # 4670.XX
Pad	VE 450 mm	VE 450 mm	VE 450 mm	VE 450 mm
Articulated Arm	6600-46	4600-46	6600-56	4600-70
Column with stop lid	-	-	Ø18 x 215 mm	Ø16 x 215 mm
Base component	Rail clamping base	Rail clamping base	Radial setting clamp	Radial setting clamp Fisso (partly optional)
Fastening straps (2 pcs)	optional Art. # 5.890-VE	optional Art. # 5.890-VE	optional Art. # 5.890-VE	optional Art. # 5.890-VE

VE = anatomically shaped visco-elastic pad

Before delivery to the clients, the components are rigidly threaded to the articulated arm and additionally fastened with glue.

Installation and use

First check

Check the product immediately after receipt for eventual transport damages and for completeness. Complaints can only be considered if the manufacturer, distributor or carrier is immediately notified through damage protocol.

Warranty

The manufacturer assumes warranty for materials and manufacturing defects as well as proper functioning of FISSO holding systems for 24 months from date of purchase.

Not covered by the warranty are normal wear, components with a shorter life time (e.g. arm pads), the consequences of inappropriate handling, combinations with external products, or damage caused by the user. Other or further claims, in particular any claims for damage are excluded. The warranty applies only if the Baitella AG repair service is used. Please see also the Sales and General Delivery Conditions.

Assembly


WARNING

All handling must be done with as little force as possible and as much force as necessary!

WARNING!

The assembly and handling of the products must be done manually and without additional tools.

Arm support lateral:

Attach the base to the side rail of the surgical table and secure it by turning the fastening screw (Pos. 1) from below clockwise. Please note that the top hook (Pos. 2) grasps behind the rail. Check if the base is fixed securely and firmly on the side rail.

Arm support vertical:

Mount the base in the same way described for arm support lateral, see above. Remove the stop lid (Pos. 3) and insert the column (Pos. 4) into the opening of the radial setting clamp. Fix it with the fastening handle (Pos. 5). Mount the stop lid. Check if the column is fixed securely and firmly into the base.

The column can be adjusted for height and inclination.

WARNING!

If the base / column is not correctly tightened, it may become loose and may cause injuries to the patient.

WARNING!


Close the gear coupling in such a manner, that the teeth grasp into each other and that they do not stand on each other.

Disassembly

The disassembly of the set must be performed also without additional tools and occurs in the reverse order of the assembly:

Hold the Arm support with one hand and loosen with the other hand the fastening screw of the rail clamping base. Remove the base from the rail of the surgical table.

Handling of the articulated arm

	1	Hold the articulated arm with one hand in the anterior segment (A) and manipulate the central clamping handle (Z) with the other hand
	2	To release turn the central handle (Z) counter clockwise (B) as far as necessary
	3	Move the articulated arm into the desired position
	4	To secure , turn the central handle (Z) clockwise (C)
	5	Check that the articulated arm is fixed firmly and is functioning properly

DANGER!

The articulated arm can move unintentionally if the central clamping is released. Always hold the arm board with one hand and manipulate the central clamping device (Z) with the other hand.

DANGER!

If the articulated arm is not tightened properly it can become loose and move, potentially causing injuries. The articulated arm must be tightened properly.

DANGER!

A damaged product can cause injuries. Use only products in proper condition and check functionality.

DANGER!

Threat through inappropriate handling
Follow instructions for use of arm support and surgical table

WARNING!

Please check the risk of collision during patient transport with exposed arm support.
If possible position the arm support parallel to the surgical table.

WARNING!

The arm support **lateral 4646.XX** has a **maximal carrying capacity of 6 kg** during intended use, the arm support **lateral 6646.XX** one of **10 kg**, the arm support **vertical 6656.XX** has a **maximal carrying capacity of 6 kg** and the arm support **vertical 4670.XX** has a **maximal carrying capacity of 4 kg**.
The maximal carrying capacity must not be exceeded.

WARNING!

The articulated arm can transmit electrical current and heat to the patient. Avoid the contact between the articulated arm and any sources of electrical current or heat.

ATTENTION!

The fixation of the articulated arm is based on the principle of friction. Changing the position without loosening the clamping mechanism can cause damage and will shorten the lifespan of the articulated arm. The articulated arm can be adjusted with little force. If the central clamping handle (Z) has been fully loosened it has to be turned in the clockwise direction!

Handling of Arm supports

The patent arm must be secured safely. Use fastening straps to secure the patient's arm safely on the arm boards.

DANGER!

If the patient's arm is not fixed properly it can come loose and move, which can cause injuries. Check secure arm positioning and fixation.

Processing (Cleaning and disinfection)

All products must be cleaned and disinfected prior to each application; this is also required for initial use after delivery (cleaning and disinfection after removal of the protective packaging). Effective cleaning and disinfection is an indispensable requirement for a safe application of the products.

You are responsible for the cleanliness of the products. Please pay attention to avoid severe contamination of the product during application (especially by handling with contaminated gloves) and follow the procedures below.

Additionally, please observe the legal provisions valid for your country as well as the hygienic instructions of the doctor's practice or of the hospital. This applies particularly to the different guidelines regarding the inactivation of prions (not relevant for the US).

Pre-treatment

Please remove coarse impurities of the products directly after use (within a maximum of 2 hrs). To this end, wipe all visually contaminated surfaces of the products intensively with disinfection wipes¹.

¹ Please consider that the wipes should be alcohol- and aldehyde-free (otherwise there is a risk of fixation of blood impurities or damage to the product), possess a fundamentally approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), be suitable for disinfecting instruments made of metallic or plastic material and compatible with the products (see chapter "Material resistance"). Please consider that a disinfectant used in the pre-treatment step is solely for the safety of the staff and cannot replace the disinfection step later to be performed after cleaning.

Manual cleaning and disinfection by wiping

Observe the following points when selecting the wipes for cleaning and disinfection:

- Fundamental suitability for the cleaning and disinfection of instruments made of metallic or plastic material
- Approved efficiency (for example VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking)
- Compatible with the products (see chapter Material resistance)

Pay attention to the instructions of the wipe manufacturer regarding application and contact time.

Proof for the general suitability of the articulated arm for an effective cleaning and disinfection was provided by an independent test laboratory using the cleaning and disinfecting mikrozid sensitive wipes (Schülke & Mayr GmbH, Norderstedt) was applied.

- Procedure:
1. The central handle of the articulated arm must be tightened during cleaning / disinfection.
 2. Wipe all surfaces of the product completely and carefully with fresh wipes. Visibly dirty or dry wipes must not be used further.
 3. Check all surfaces for remaining visible contamination and wipe again by applying a fresh wipe, if required.
 4. Check the products (see chapters "Control" and "Maintenance").
 5. Reassemble the product and wipe again using a fresh wipe.

DANGER!

The product will be delivered non sterile and is **non-sterilizable**. Prior to the initial and each following use, the product has to be cleaned and disinfected as well as checked for visual irregularities and malfunctioning according to the indications given in this manual.

ATTENTION!

Material damage due to improper cleaning and disinfection.

This product can be cleaned and disinfected **only manually**. The cleaning and disinfection agents must be compatible with the product.

Metallic surface:

Use only cleaning and disinfection wipes. Never use metal brushes or steel wool.

Integral foam pad made of polyurethane:

Surface disinfectants containing a combination of the active ingredients quaternary compounds and guanidine derivatives are suitable.

Visco-elastic pads

Rinse with clean water and dry with soft fabric.

The following products are strictly forbidden:

- undiluted alcohol and / or acetone
- Perchlorethylene, Trichlorethylene and wax
- abrasive cleaning products

Mechanically damaged pads with tears and cuts have to be replaced with new pads to avoid the entry of moisture (hygienic requirements no longer fulfilled).

WARNING!

Danger of explosion.

Alcohol containing agents form ignitable mixtures, which can explode if used together with high-frequency applications.

Do not use alcohol containing agents together with high-frequency applications.

WARNING!

The articulated arm must not be sprayed directly with liquid or immersed in liquids.

The central handle of the articulated arm must be tightened during cleaning / disinfection.

Control

Check after cleaning/disinfection the action of moving parts to ensure smooth operation throughout the intended range of motion, the clamping force, the surface regarding corrosion, damages, chipping, and dirt and sort out impaired products. Products that are still dirty must be cleaned and disinfected again.

Packaging and storage

WARNING!

Store the articulated arms such that they are not damaged.

Please avoid re-contamination of the product by adequate protection / packaging and store the products in a dry and dust-free place. After assembly and prior to next use, further disinfection by wiping all surfaces (as in pre-treatment) may be required.

Maintenance

Special maintenance is not required. The product should be sent back to the manufacturer or distributor if the action of moving parts or the clamping force is impaired, or in case of damage.

Instrument oil may never be used.

NOTE

Repair must be carried out by Baitella AG. If not respected, the warranty becomes void!

ATTENTION!

Please send back only reprocessed products to the distributor / manufacturer for repair (cleaned and disinfected).

Material resistance

Please note that the substances listed are **not** ingredients of the cleaning or disinfection agent:

- Organic, inorganic and oxidizing acids (permitted pH higher than 6.5)
- Bases (permitted pH lower than 8.5, neutral detergents are recommended)
- Organic compounds / solvents (for example: alcohol, aldehyde, acetone, ether, benzene)
- Oxidizing agents (for example: peroxide)
- Halogens (chlorine, iodine, bromine)
- Aromatic, halogenated hydrocarbons
- Wax

The product is non-sterilizable. Please do not expose any products to temperatures higher than 60 °C (140 °F).

Reusability

The products can be reused if processed properly and undamaged (see chapter "Control"). The user assumes full responsibility if damaged or dirty products are used (no liability will be accepted if this is not observed).

Obligation to Notify

The purchaser shall inform the manufacturer about risks and recalls, as well as about reportable incidents, insofar as a product supplied by manufacturer is affected.

Disposal

The products must be correctly disposed of in accordance with national regulations and medical guidelines.

Technical data

	Arm support lateral		Arm support vertical	
	4646.XX	6646.XX	4670.XX	6656.XX
Classification according to Medical Device Directive (EU) 2017/745	class 1	class 1	class 1	class 1
Classification FDA	class 1	class 1	class 1	class 1
Radius of movement of the articulated arm	460 mm	460 mm	700 mm	560 mm
Column length	-	-	215 mm	215 mm
Total length	approx. 580 mm	approx. 580 mm	approx. 1030 mm	approx. 877 mm
Maximal carrying capacity of the articulated arm during intended use	6 kg	10 kg	4 kg	6 kg
Weight	approx. 3.9 kg according to model	approx. 4.3 kg according to model	approx. 4.1 kg according to model	approx. 5.4 kg according to model

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