

**Simply the best fixation.**



**FISSO**  
Swiss Made

## Patient positioning / arm support system



### Features:

- « one step set up »
- Central clamping mechanism for easy, fast and secure positioning of the arm support
- Rapid locking base allows connection to the side rails of the OP table
- Uncompromising, efficient handling and safe holding
- Space saving on the OR table, compact design
- Anatomically shaped visco-elastic (VE) – pad



**FISSO**  
Swiss Made

**Simply the best fixation.**

## Patient arm support lateral



**6646.01**

| Order no.      | Description                                  | Pad           | Articulated arm |        | Base Order no. | Weight  |
|----------------|--|---------------|-----------------|--------|----------------|---------|
|                |  |               | Order no.       | Length |                |         |
| <b>6646.01</b> | Arm support lateral<br>6600-46, strong model | visco-elastic | 2.270<br>460 mm | 10 kg  | 5.751          | 4.30 kg |

Subject to modifications

## Patient arm support vertical



**6656.02**

| Order no.      | Description                                   | Pad           | Articulated arm |        | Base Order no. | Weight  |
|----------------|---|---------------|-----------------|--------|----------------|---------|
|                |   |               | Order no.       | Length |                |         |
| <b>6656.02</b> | Arm support vertical<br>6600-56, strong model | visco-elastic | 2.275<br>560 mm | 6 kg   | 5.750          | 5.40 kg |

Subject to modifications

## Additional components



**5.890-VE**

| Order no.       | Description      | Length | Width | Quantity | Weight  |
|-----------------|------------------|--------|-------|----------|---------|
| <b>5.890-VE</b> | Fastening straps | 600 mm | 48 mm | 2 pcs    | 0.06 kg |

Subject to modifications

Baitella AG  
Thurgauerstrasse 70, CH-8050 Zurich, Switzerland  
Phone +41 44 305 80 00, Fax +41 44 305 80 05  
info@baitella.com, www.fisso.com



BMPASe 02/20